

APC Anti-Human CD27 (O323)

Catalog Number: APC-65132

Basic Information

Catalog Number:

APC-65132

Size:

100 tests, 5 µl / test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC012160

GeneID (NCBI):

939

ENSEMBL Gene ID:

ENSG00000139193

UNIPROT ID:

P26842

Full Name:

CD27 molecule

Calculated MW:

29 kDa

Purification Method:

Affinity purification

CloneNo.:

O323

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

CD27 (also known as TNFRSF7) is a type I glycoprotein expressed on some B cells and the majority of T cells, and is a member of the tumor necrosis factor (TNF) receptor family. CD27 is required for generation and long-term maintenance of T cell immunity (PMID: 11062504). It is a receptor for CD70 (CD27L). Ligation of CD27 by CD70 induces strong ubiquitination of TRAF and the activation of both canonical and non-canonical NF-kappaB pathways, as well as the JNK pathway (PMID: 19426224). CD27 may also play a role in apoptosis through association with SIVA1.

Storage


Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.


Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 5 ul APC Anti-Human CD27 (O323) (APC-65132, Clone:O323) or APC Mouse IgG1 Isotype Control (MOPC-21) (APC-65124, Clone: MOPC-21), and 5 ul CoraLite® Plus 488 Anti-Human CD19 (HIB19) (CL488-65110, Clone: HIB19). Cells were not fixed.



1x10⁶ human PBMCs were surface stained with 5 ul APC Anti-Human CD27 (O323) (APC-65132, Clone:O323) or APC Mouse IgG1 Isotype Control (MOPC-21) (APC-65124, Clone: MOPC-21). Cells were not fixed.