

For Research Use Only

APC Anti-Human CD27 (O323)

Catalog Number: APC-65132



Basic Information

Catalog Number: APC-65132	GenBank Accession Number: BC012160	Purification Method: Affinity purification
Size: 100tests, 5 µl/test	GeneID (NCBI): 939	CloneNo.: O323
Source: Mouse	ENSEMBL Gene ID: ENSG00000139193	Excitation/Emission maxima wavelengths: 650 nm / 660 nm
Isotype: IgG1, kappa	UNIPROT ID: P26842	
	Full Name: CD27 molecule	
	Calculated MW: 29 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

CD27 (also known as TNFRSF7) is a type I glycoprotein expressed on some B cells and the majority of T cells, and is a member of the tumor necrosis factor (TNF) receptor family. CD27 is required for generation and long-term maintenance of T cell immunity (PMID: 11062504). It is a receptor for CD70 (CD27L). Ligation of CD27 by CD70 induces strong ubiquitination of TRAF and the activation of both canonical and non-canonical NF-kappaB pathways, as well as the JNK pathway (PMID: 19426224). CD27 may also play a role in apoptosis through association with SIVA1.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human peripheral blood lymphocytes were surface stained with 5.00 ul/test APC Anti-Human CD27 (APC-65132, Clone:O323) (red) or isotype control antibody (blue). Cells were not fixed.