

For Research Use Only

APC Anti-Human CD137 (4B4-1)

Catalog Number: **APC-65159**



Basic Information

Catalog Number:

APC-65159

Size:

100 tests, 10 µl/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC006196

GeneID (NCBI):

3604

UNIPROT ID:

Q07011

Full Name:

tumor necrosis factor receptor
superfamily, member 9

Calculated MW:

28 kDa

Purification Method:

The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

4B4-1

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

Human, Non-human primates

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ PHA stimulated (5 μ g/mL, 3 days) (red) or unstimulated (blue) human peripheral blood mononuclear cells (PBMCs) were surface stained with 10.00 μ l APC Anti-Human CD137 (APC-65159, Clone: 4B4-1). Cells were not fixed.