

For Research Use Only

APC Rat IgG2a Isotype Control (2A3)

Catalog Number: APC-65209



Basic Information

Catalog Number:

APC-65209

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Affinity purification

Size:

100ug, 0.2 mg/mL

Full Name:

CloneNo.:

2A3

Source:

Rat

Excitation/Emission maxima
wavelengths:

650 nm / 660 nm

Isotype:

IgG2a, kappa

Applications

Tested Applications:

FC

Species Specificity:

n/a

Background Information

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.1 ug APC Rat IgG2a Isotype Control (2A3) (APC-65209, Clone:2A3) (red), or 0.1 ug APC Anti-Mouse CD4 (RM4-5) (APC-65141) (green), or unstained (blue). Cells were not fixed.