

For Research Use Only

APC Anti-Human SEMA7A Rabbit Recombinant Antibody

Catalog Number: APC-98324-3



Basic Information

Catalog Number:

APC-98324-3

Concentration:

100 tests, 5 µl/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2005

GenBank Accession Number:

NM_003612.5

GeneID (NCBI):

8482

UNIPROT ID:

O75326-1

Full Name:

semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)

Calculated MW:

75kDa

Purification Method:

Protein A purification

CloneNo.:

241814B12

Recommended Dilutions:

FC: 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : PHA treated human PBMCs,

Background Information

SEMA7A, also named CD108 and SEMA7, belongs to the semaphorin family. SEMA7A is the only membrane-associated glycosylphosphatidylinositol (GPI)-linked semaphorin and plays a role in the regulation of neuronal axon outgrowth. In addition to its role in guiding axon pathfinding during neuronal development, SEMA7A has diverse functions in morphogenesis and immune cell control and regulates T cell responses via the $\alpha 1 \beta 1$ integrin receptor.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or PHA-treated human PBMCs were surface stained with 5 ul APC Anti-Human SEMA7A Rabbit RecAb (APC-98324-3, Clone: 241814B12). Cells were not fixed.



1x10⁶ PHA-treated human PBMCs were surface stained with 5 ul APC Anti-Human SEMA7A Rabbit RecAb (APC-98324-3, Clone: 241814B12) (red) or APC Rabbit IgG Isotype Control Recombinant Antibody (APC-98136, Clone: 240953C9) (blue). 1x10⁶ untreated human PBMCs were surface stained with 5 ul APC Anti-Human SEMA7A Rabbit RecAb (APC-98324-3, Clone: 241814B12) (black, dashed). Cells were not fixed.