

For Research Use Only

FcZero-rAb™ APC Anti-Human TLR3/CD283 Rabbit Recombinant Antibody

Catalog Number: APC-FcA98016



Basic Information

Catalog Number:

APC-FcA98016

Concentration:

100 tests, 5 µl/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34044

GenBank Accession Number:

BC094737

GeneID (NCBI):

7098

Full Name:

toll-like receptor 3

Calculated MW:

904 aa, 104 kDa

Purification Method:

Protein A purification

CloneNo.:

230445F1

Recommended Dilutions:

FC (Intra): 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Positive Controls:

FC (Intra): human immature monocyte-derived dendritic cells,

Background Information

TLR3 belongs to the Toll-like receptor family which is important in the innate immune response to pathogens. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. TLR3 is a nucleotide-sensing TLR that is activated by double-stranded RNA. Upon recognition of dsRNA, TLR3 transmits signals via the adaptor protein TICAM-1/TRIF, leading to the induction of type I IFN, cytokine/chemokine production, and dendritic cell (DC) maturation (PMID: 18262679).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human immature monocyte-derived dendritic cells were intracellularly stained with 5 μ L APC Anti-Human TLR3/CD283 Rabbit RecAb (APC-FcA98016, Clone: 230445F1) or APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9). Cells were co-stained with Anti-Human CD11c. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).