

For Research Use Only

# FcZero-rAb® APC Anti-Human CDCP1/CD318 Rabbit Recombinant Antibody

Catalog Number: APC-FcA98294



## Basic Information

Catalog Number:

APC-FcA98294

Concentration:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2574

GenBank Accession Number:

BC153169

GeneID (NCBI):

64866

UNIPROT ID:

Q9H5V8

Full Name:

CUB domain containing protein 1

Purification Method:

Protein A purification

CloneNo.:

242428F9

Recommended Dilutions:

FC: 5 ul per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

## Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : HT-29 cells,

## Background Information

CDCP1, also called CD318, TRASK, and SIMA 135, is a 135-140 kDa transmembrane protein containing three CUB domains within the extracellular region and a hexalysine stretch within the cytoplasmic region. A soluble form with MW 80 kDa (p80) is cleaved at an unidentified site on the N-terminal side. CDCP1 is a novel marker for hematopoietic stem cells.

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> HT-29 cells (left) or MCF-7 cells (right) were stained with 5 ul APC Anti-Human CDCP1/CD318 Rabbit RecAb (APC-FcA98294, Clone: 242428F9) (red) or APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.