

For Research Use Only.
Not For Use In Diagnostics.

APC-Cyanine7-conjugated Streptavidin Rabbit RecAb

Catalog Number: AY7-PF00030



Information

Catalog Number:

AY7-PF00030

Host:

Applications:

FC, IF

Reactivity:

Physical State:

Liquid

Conjugation:

APC-Cyanine7

Recommended Dilutions

1:200-1:1000 for flow cytometry and immunofluorescence

Storage

Storage:

2-8°C for one year.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.4.

Aliquoting is unnecessary

Purity

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were unstained (blue) or surface stained with Biotin Anti-Human CD8 (UCHT4) Mouse IgG2a Recombinant Antibody (Biotin-65559, Clone: UCHT4) (red), washed with PBS, and stained with Streptavidin-CorLite® Plus 405 Conjugate (Cat.No. AY7-PF00030) and analyzed by flow cytometry. Cells were not fixed. Lymphocytes were gated.