

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ITGA6 Fusion Protein



## Basic Information

**Catalog Number:**

Ag0527

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

FNLDTREDNVIKRYGDPGSLFGFSLAMHWQLQPEDKR  
LLLVGAPRAEALPLQRANRTGGLYSCDITARGPCTRIE  
FDNDADPTSESKEDQWMGVTVQSQGGGKVVTC  
HRYEKRQHVNTKQESRDIFGRCYVLSQNLRIEDDMDG  
GDWSFCDGRLRGHEKFGSCQQGVAATFTKDFHYIVF  
GAPGTYNWKGIVRVEQKNNTFFDMNIFEDGPYEVGG  
ETEHDSESLVPVPANSYLGFLDSGKGVSKDEITFVSG  
APRANHSGAVVLLKRDMSAHLLEPHIFDGEGLASSF  
GYDVAVVLDLNKDGWQDIVIGAPQYFDRDGEVGGAV  
YYMNNQQGRWNNVKPIRLNGTKDSMFGIAVKNIGDI  
NQDGYPDIAVGAPYDDLKGVFIYHGSANGINTKPTQ  
VLKGISPYFGYSIAGNMDLDRNSYPDVAVGSLSDSVTI  
FRSRPVINIQTITVTPNRI DLRQKTACGAPSGICLQV  
KSCFEYTANPAGYNPSISIVGTLE  
(24 - 523 aa encoded by NM\_001079818)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.