

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Sf3b3 Fusion Protein



Basic Information

Catalog Number:

Ag0544

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

KVHTLLTVEVFQVIRSLMAFRLTGGTKDYIVVGSDSGR
IVILEYQPSKNMFEKIHQETFGKSGCRRIVPGQFLAVD
PKGRAVMISAIKQKLVYILNRDAAARLTSSPLEAHKA
NTLVYHVVGVDVGFENPMFACLEMDYEEADNDPTG
EAAANTQQLTFYELDLGLNHVVRKYSEPLEEHGNFLI
TVPGGSDGPGSVLICSENYITYKNFGDQPDIRCPIRR
RNDLDDPERGMIFVCSATHKTKSMFFFLAQTEQGDIF
KITLETDEDMVTEIRLKYFDTVPVAAAMCVLKTGFLEV
ASEFGNHLYQIAHLGDDDEEPEFSSAMPLEEGDTFF
FQPRPLKNLVLVDELDSLPLFCQIADLANEDTPQLY
VACGRGPRSSLRVLRHGLEVSEMAVSELPGNPNVAVW
TVRRHIEDEFDAYIIVSFVNATLVLSIGETVEEVTDSGF
LGTTPTLSCSLLGDDALVQVYPDGI
(49-525 aa encoded by BC139016)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.