

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Integrin alpha-5 Fusion Protein



Basic Information

Catalog Number:

Ag0860

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

FTVPHLRDTKKTIQDFQILSKNLNNSQSDVVVSRFSV
EAQAQVTLNGVSKPEAVLFPVSDWHPRDQPQKEED
LGPAVHHVYELINQGPSSISQGVLELSCPQALEGQQL
LYVTRVTGLNCTTNHPINPKGLELDPEGSLHHQKRE
APSRSSASSGPQILKCPAEAEFRLRCELGPLHQEQESQ
SLQLHFRVWAKTFLQREHQPFSLQCEAVYKALKMPY
RILPRQLPQKERQVATAVQWTKAEGSYGVPLWIIILAI
LFGLLLLGLLIYLYKLGFFKRSLPYGTAMEKAQLKPPAT
SDA

(749-1049 aa encoded by BC008786)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.