

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant mouse Ttll5 protein



Basic Information

Catalog Number:

Ag0919

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

e coli-derived, PKG, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MARDLEETASSEDEDLANQEDHPCIMWTGGCRRIP
VLVFAEAAILTKDNNIRVIGERYHLSYKIVRTDSRLVRSI
LTAHGFHEVHPSSTDYINLMWTGSHLKPFLLRTLSEAQ
KVNHFPRSYELTRKDRLYKNIIRMQHTHGFKAFHILPQ
TFLLPAEYAEFCNSYSKDRGPWIVKPVASSRGRGVYLI
NNPNQISLEENILVSRVINNPLIIDDFKFDVRLYVVLVTS
YDPLVIYLYEEGLARFATVRYDQGSKNIRNQFMHLTNY
SVNKKSGDYVSCDDPEVEDYGNKWSMSAMLRYLKQE
GKDTTALMAHVEDLIIKTIISAELAIATACTFVPHRSS
CFELYGFDVLIIDNTLKPWLLVNLSPSLACDAPLDLKK
ASMSDMFTVVGVFCQDPAQ
(1-400aa encoded by BC035276)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.