

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FOXP4 protein



Basic Information

Catalog Number:

Ag10034

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

NTEHALDDRSTAQCRVQM~~QVVQ~~QLEIQLAKESERL
QAMMAHLHMRPSEPKPFSQPLNPVPGSSSF~~SKVTVS~~
AADSFPDGLVHPPTSAAAAPVTPLRPPGLGASLHG~~GG~~
PARRRSSDKFCSPISSELAQNHEFYKNADV~~RPPFTYAS~~
LIRQAILET~~PDRQLTLNEIYNWFTRMFAYFRRNTATWK~~
NAVRHNLSLHKCFVRVENVKGAVVTV~~DEREYQKRRP~~
PKMTGSPTLVKNMISGLSYGALNASYQAALA~~ESSFLL~~
NSPGMLNPGSASSLLPLSHDDVGAPVE~~PLPSNGSSSP~~
PRLSPPQYSHQVQVKEEPAEA~~EEDRQPGPLGAPNP~~
SASGPPEDRDLEELPGEELS

(327-678 aa encoded by BC052803)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.