

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# GALNT14 Fusion Protein



## Basic Information

Catalog Number:

Ag10344

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e. coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EVNRDWLQPLLHRVKEDYTRVVCVIDIINLDTFTYIE  
SASELRGGFDWSLHFQWEQLSPEQKARRLDPTPIRT  
PIIAGGLFVIDKAWFDYLGKYDMDMDIWGGENFEISF  
RVWMCSSLEIVPCSRVGHVFRKKHPYVFPDGNANT  
YIKNTKRTAEVWMDEYKQYYAARPFALERPFGNVES  
RLDLRKNLRCQSFKWYLENIPELSIPKESSIQKGNIRQ  
RQKCLEQRQNNQETPNLKLSPCAKVKGEDAKSQV  
WAFYTYQILQEELCSVITLFPGAPVVLVLCCKNGDD  
RQQWTKTGSHEIHASHLCLDMDMFGDGTENGKEIV  
VNPCESSLMSQHWDMVSS

(203-552aa encoded by BC010659)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.