

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

EED Fusion Protein



Basic Information

Catalog Number:

Ag10358

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSEREVSTAPAGTDMPAAKKQKLSSDSENSNPDLSGD
ENDDAVSIESGTINTERPDTPTNPNAPGRKSWGKGGK
WKSkkCKYSFKCVNSLKEDHNQPLFGVQFNWHSKEG
DPLVFATVGSNRVTLYECHSQGEIRLLQSYVDADADE
NFYTCAWTYDSNTSHPLLA VAGSRGIIIRINPITMQCI
KHVYVGHGNAINELKFHPRDPNLLSVSKDHALRLWNI
QTDTLVAIFGGVEGHRDEVL SADYDLLGEKIMSCGMD
HSLKLWRIN SKRMMNAIKESYDYNPNKTRPFISQKIH
FPDFSTRDIHRNYVDCVRWLGDLILSKSGRAILHSHQ
QCMRDPVSPNLRQHLSCE NAIVCWKPGKMEDDIDKI
KPSESNVTILGRFDYSQCDIWYMRFSMDFWQKMLAL
GNQVGKLYVWDEVEDPHKAKCTTLTHHKCGAAIRQ
TSFSRDSIILJAVCDDASIRWRDLR
(1-466 aa encoded by BC068995)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.