

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

ACAD10 Fusion Protein



Basic Information

Catalog Number:

Ag10529

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MGGVLI P S P G R V A A E W E V Q N R I P S G T I L K A L M E G G E N
G P W M R F M R A E I T A E G F L R E F G R L C S E M L K T S V P V D S F
F S L L T S E R V A K Q F P V M T E A I T Q I R A K G L Q T A V L S N N F Y
L P N Q K S F L P L D R K Q F D V I V E S C M E G I C K P D P R I Y K L C L E
Q L G L Q P S E S I F L D D L G T N L K E A A R L G I H T I K V N D P E T A
V K E L E A L L G F T L R V G V P N T R P V K K T M E I P K D S L Q K Y L K
D L L G I Q T T G L

(1-237 aa encoded by BC015056)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.