

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PTPN18 protein



## Basic Information

Catalog Number:

Ag10610

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSRSLDSARSFLERLEARGGREGAVLAGEFSKRCERY  
WAQEQEPLQTGLFCITLIKEKWLNEDIMLRCLKVTFQK  
ESRSVYQLQYMSWPDGRGVPSSPDHMLAMVEEARRL  
QSGGPEPLCVHCSAGCGRTGVLCVLDYVRRQLLTTQM  
IPPDFSLFDVVLKMRKQRPAAVQTEEQYRFLYHTVAQ  
MFCSTLQNASPHYQNIKENCAPLYDDALFLRTPQALL  
AIPRPPGGVLRISIVPGSPGHAMADTYAVVQKRGAP  
AGAGSGTQTGTGARSAAEAPLYSKVTPRAQRPGA  
HAEDARGTLPGRVPADQSPAGSGAYEDVAGGAQTG  
GLGFNLIRGPKGRDPPAEWTRV  
(1-351 aa encoded by BC052800)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.