

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human RSNB1 protein



## Basic Information

**Catalog Number:**

Ag10628

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MAAQVGA VRVRAVAAQEEPDKEGKPKPHAGVSPR  
GVKRRSSSSGGSQEKGRPSQEPPLAPPHRRRSRQ  
HPGPLPPTNAAPTVPVPEPLPPPPPSLAPAGPAV  
AAPLPAPSTSALFTFSPLTVSAAGPKHKHGERHKKH  
HHRGPDGDPSSCGTDLKHKDKQENGERTGGVPLIKAP  
KRETPDENGKTQRADDFVLKKIKKKKKKHREDMRGR  
RLKMYNKEVQTV CAGLTRISKEILTQGGQINSTSGLNKE  
SFRYLKDEQLCRLNLGMQEYRV PQGVQTPFMTHQEH  
SIRRNFLKTGTFKSNFIHEEHQSNGGALV LHAYMDELS  
FLSPMEMERFSEEFALTFSENEKNAAYALAI VHGAA  
AYLPDFLDYFAFNFPNTPVKMEILGKKDIETTTISNFHT  
Q

(1-411aa encoded by BC026155)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.