

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

DTNA Fusion Protein



Basic Information

Catalog Number:

Ag1075

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MIEDSGKRGNTMAERRQLFAEMRAQDLDRIRLSTYRT
ACKLRFVQKICNLHLVDIWNVIEALRENALNNDPNT
ELNVSRLA V LSTIFYQLNKRMPPTHQIHVEQSISLLN
FLAAFDPEGHGKISVFAVKMALATLCGGKIMDKLRYI
FSMISDSSGVMVYGRYDQFLREV LKLP TAVFEGPSFG
YTEQSARSCFSQQKKT V L N G F L D T L M S D P P P Q C L V W
L P L L H R L A N V E N V F H P V E C S Y C H S E S M M G F R Y R C Q Q
C H N Y Q L C Q D C F W R G H A G G S H S N Q H Q M K E Y T S W K S
P A K K L T N A L S K S L S C A S S R E P L H P M F P D Q P E K P L N L A H
I V P P R P V T S M N D T L F S H S V P S S G S P F I T R S S D G A F G G C
V

(1-371 aa encoded by BC005300)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.