

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PDHA2 Fusion Protein



Basic Information

Catalog Number:

Ag10882

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MLAAFISRVLRRVAQKSARRVLVASRNSNDATFEIKK
CDLYLLEEGPPVTTVLTRAEGLYRMMMLTVRRMELKA
DQLYKQKFIRGFCHLCDGQEACCVGLEAGINPSDHVI
TSYRAHGVCYTRGLSVRSILAEELTGRGGCAKGGGS
MHMYTKNFYGGNGIVGAQGPLGAGIALACKYKGNDE
ICLTLYGDGAANQGQIAEAFNMAALWKLPCVFICEN
NLYGMGTSTERAAASPDYYKRGNFIPGLKVDGMDVLC
VREATKFAANYCRSGKGPILMELQTYRYHGHSMSPDG
VSYRTREEIQEVRSKRDPIILQDRMVNSKLATVEELKE
IGAEVRKEIDDAQFATTDPEPHLEELGHIIYSSDSSF
EVRGANPWIKFKSVS

(1-388 aa encoded by BC030697)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.