

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HLA class I (HLA-B) Fusion Protein



## Basic Information

**Catalog Number:**

Ag10947

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MRVTAPRTVLLLLSGALALTETWAGSHSMRYFYTAMS  
RPGRGEPRFISVGYVDDTQFVRFSDAASPREEPRAP  
WIEQEGPEYWRNTQICKTNTQTYRESLRNLRGYYN  
QSEAGSHTLQRMYGCDVGPDRLLRGHDQYAYDGK  
DYIALNEDLSSWTAADTAAQITQRKWEAAREAEQLR  
AYLEGLCVEWLRRYLENGKETLQRADPPKTHVTHHPIS  
DHEATLRCWALGFYPAEITLTWQRDGEDQTQDTELV  
ETRPAGDRTFQKWAAVVVPSPGEEQRYTCHVQHEGL  
PKPLTLRWEPSQSTIPIVGIVAGLAVLAVVVIGAVVA  
TVMCRRKSSGGKGGSYSQAASSDSAQGSQSDVSLTA  
(1-362 aa encoded by BC013187)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.