

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KLHDC2 protein



## Basic Information

**Catalog Number:**

Ag10951

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MADGNEDLRADDLPGPAFESYESMELACPAERSGHV  
AVSDGRHMFVWGGYKSNQVRGLYDFYLPREELWIYN  
METGRWKKINTEGDVPPMSGSCAVCVDRVLYLFGG  
HHSRGNTNKFYMLDSRSTDRVLQWERIDCQGIPPSSK  
DKLGVVWYKNKLIFGGYGYLPEDKVLGTFEFDTSF  
WNSSHPRGWNDHVHILDTETFTWSQPITTKAPSPR  
AAHACATVGNRGFVFGGRYRDARMNDLHYLNLDTW  
EWNELIPQGI CPVGRSWHSLTPVSSDHLFLFGGFTTD  
KQPLSDAWTYCISKNEWIQFNHPYTEKPRLWHTACA  
SDEGEVIVFGGCANNLLVHHRAAHSNEILF SVQPKSL  
VRLSLEAVICFKEMLANSWNCLPKHLLHSVNQRFGSN  
NTSGS

(1-406 aa encoded by BC024192)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.