

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human HMGCCL1 protein



## Basic Information

### Catalog Number:

Ag10985

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

MGNVPSAVKHCLSYQQLLREHLWIGDSVAGALDPAQ  
ETSQLSGLPEFVKIVEVGPRDGLQNEKVIIVPTDIKIEFI  
NRLSQTGLSVIEVTSFVSSRWVPQADHTEVMKGIH  
QYPGVRYPVLTPNLQGFHHAAGATEISVFGAASE  
SFSKKNINCSEESMGKFEEVVKARHMNIPARGYVS  
CALGCPYEGSITPQKVTEVSKRLYGMGCYEISLGDITG  
VGTPGSMKRMLESVMKEIPPGALAVHCHDTYGQALA  
NLTALQMGINVVDSAVSGLGGCPYAKGASGNVATE  
DLIYMLNGLGLNTGVNLYKMEAGDFICKAVNKTNS  
KVAQASFNA

(1-340 aa encoded by BC024194)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.