

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RB1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag11035

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

SKDREGPTDHLSEACPLNPLQNNHTAADMYLSPVRS  
PKKKGSTTRVNSTANAETQATSAFQTQKPLKSTLSL  
FYKKVYRLAYLRNLTCERLLSEHPELEHIIWTLFQHTL  
QNEYELMRDRHLDQIMMCSMYGICKVKNIDLKFKIIVT  
AYKDLPHAVQETFKRVLIKEEYDSIIVFYNSVFMQRL  
KTNILQYASTRPPTLSPIPHIPRSPYKFPSSPLRIPGGNIY  
ISPLKSPYKISEGLPTPTKMTPRSRI LVSIGESFGTSEKF  
QKINQMVCNSDRVLKRSAEGSNPPKPLKLRFDIEGS  
DEADGSKHLPGESKQQKLAEMTSTRTRMQKQKMN  
DSMDTSNKEEK

(576-928 aa encoded by BC039060)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.