

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human KRT16 protein



Basic Information

Catalog Number:

Ag11240

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MQNLNDRLASYLDKVRAL E E A N A D L E V K I R D W Y Q R Q
R P S E I K D Y S P Y F K T I E D L R N K I I A A T I E N A Q P I L Q I D N A R
L A A D D F R T K Y E H E L A L R Q T V E A D V N G L R R V L D E L T L A R
T D L E M Q I E G L K E E L A Y L R K N H E E M L A L R G Q T G G D V N
V E M D A A P G V D L S R I L N E M R D Q Y E Q M A E K N R R D A E T
W F L S K T E E L N K E V A S N S E L V Q S S R S E V T E L R R V L Q G L E
I E L Q S Q L S M K A S L E N S L E E T K G R Y C M Q L S Q I Q G L I G S V
E E Q L A Q L R C E M E Q S Q S Q E Y Q I L L D V K T R L E Q E I A T Y R R
L L E G E D A H L S S Q Q A S G Q S Y S S R E V F T S S S S S S R Q T R
P I L K E Q S S S S F S Q G Q S S

(121-473 aa encoded by BC039169)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.