

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# C1R Fusion Protein



## Basic Information

**Catalog Number:**

Ag11242

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

AVCQDDGTWHRAMPRCKIKDCGQPRNLPNGDFRYT  
TTMGVNTYKARIQYYCHEPYKMQTRAGSRESEQGV  
YTCTAQGIWKNEQKGEKIPRCLPVCCKPVNPVEQRQ  
RIIGGQKAKMGNFPWQVFTNIHGRGGGALLGDRWIL  
TAAHTLYPKEHEAQSNASLDVFLGHTNVEELMKLGNH  
PIRRVSVHPDYRQDESYNFECDIALLENSVTLGPNLL  
PICLPDNDTFYDLGLMGYVSGFGVMEEKIAHDLRFVRL  
PVANPQACENWLRGKNRMDVFSQNMFCAGHPSLK  
QDACQGDSSGGVFAVRDPNTDRWVATGIVSWGIGC  
SRGYGFYTKVLLNYVDWIKKEMEEED  
(356-705 aa encoded by BC035220)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.