

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human VPS4B protein



Basic Information

Catalog Number:
Ag11583

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MSSTSPNLQK AIDLASKAAQEDKAGNYEALQLYQH
AVQYFLHVVKYEAQGDKAKQSIRAKCTEYLDRAEKLK
EYLNKEKKAQKPVKEGQPSPADEKGNDSGEGESD
DPEKKLQNLQGAIVIERPNVKWSDVAGLEGAKEAL
KEAVILPIKPHLFTGKRTPWRGILLFGPPGTGKSYLAK
AVATEANNSTFFSISSSDLVSKWLGESEKLVKNLFQLA
RENKPSIIFIDEIDSLCGSRSENESEAARRIKTEFLVQM
QGVGVNDGILVLGATNIPWVLD SAIRRRFEKRIYIPL
PEPHARAAMFKLHLGTTQNSLTEADRELGRKTDGYS
GADISIIVRDALMQPVRKVSATHFKKVRGSPRADPN
HLVDDLLTPCSPGDPGAIEMTWMDVPGDKLLEPVVS
MSDMLRSLSNKPTVNEHDLKLLKFFTEDFGQEG
(1-444 aa encoded by BC039574)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.