

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SLC25A44 protein



## Basic Information

**Catalog Number:**

Ag11595

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MEDKRNQIIWEHLDKKKFYVFGVAMTMMIRVSVYP  
FTLIRTRLQVQK GKSLYHGTFDAFIKLRADGITGLYRG  
FLVNTFTLISGQCYVTTYELTRKRVADYSQSNVKSIV  
AGGSASLVAQSITVPI DVVSQHLMMQRKGEKMGRF  
QVRGNPEGQGVVAFGQTKDII RQILQADGLRGFYRG  
YVASLLTYIPNSAVWWPFYHFYAEQLSYLCPKECPHI  
VFQAVSGPLAAATASILTNPMDVIRTRVQVEGKNSIIL  
TFRQLMAEEGPWGLMKGLSARIISATPSTIVIVVGYES  
LKKLSLRPELVDSRHW

(1-314 aa encoded by BC039854)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.