

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human SCPEP1 protein



Basic Information

Catalog Number:

Ag11931

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KTTWLQAASLLFVDNPVGTGFSYVNGSGAYAKDLAM
VASDMMVLLKTFFSCHKEFQTVPFYIFSESYGGKMAA
GIGLELYKAIQRGTIKCNFAGVALGDSWISPVDSVLS
WGPPLYSMSLLEDKGLAEVSKVAEQVLNAVNGLYRE
ATELWGAEMII EQNTDGVNFYNI LTKSTPTSTMESSL
EFTQSHLVCLCQRHVRHLQRDALSQLMNGPIRKLKII
PEDQSWGGQATNVFVNMEEDFMKPVISIVDELLEAG
INVTVYNGQLDLIVDTMGQEAWVRKLKWPPELPKFSQ
LKWKALYSDPKSLETSAFVKS YKNLAFYWILKAGHMV
PSDQGD MALKMMRLVTQQE

(102-452 aa encoded by BC072405)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.