

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PSAP Fusion Protein



## Basic Information

Catalog Number:

Ag1194

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MYALFLLASLLGAALAGPVLGLKECTRGSAVWCQNVK  
TASDCGAVKHCLQTVWNKPTVKSLPCDICKDVVTAA  
GDMLKDNATEEEILVYLEKTCDWLPKPNMSASCKEIV  
DSYLPVILDIIKGEMSRPGEVCSALNLCESLQKHLAELN  
HQQKQLESNKIPELDMTEVVAPFMANIPLLLYPQDGPR  
SKPQPKDNGDVCQDCIQMVTDIQTAVRTNSTFVQA  
LVEHVKEECDRLGPGMADICKNYISQYSEIAIQMMMH  
MQPKEICALVGFCDEVKEMPMQTLVPAKVASKNVIP  
ALELVEPIKKHEVPAKSDVYCEVCEF  
(1-320 aa encoded by BC001503)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.