

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PPRC1 Fusion Protein



## Basic Information

**Catalog Number:**  
Ag11940

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
VIPEVGSRWNVKRHQDITIKPVLISLGPAAAPPPCIAAS  
REPLDHRTSSEQADPSAPCLAPSSLLSPEASPCRNDM  
NTRTPPEPSAKQQRSMRCYRKACRSASPSSQGWQGRR  
GRNSRSVSSGNSRTSEASSSSSSSSSSSSRSRSLSP  
HKRWRRSSCSSSGRSRRCSSSSSSSSSSSSSSSSSS  
RSRSRSPSPRRRSDRRRSYRSHDHYQRQRLQKERAI  
EERRVVFIGKIPGRMTRSELKQRFVSVFGEIECTIHFVRV  
QGDNYGFVTYRYAEEFAAIESGHKLRQADEQPFDL  
CFGRRQFCRKRSYDLDNSNREDFDPAPVKSKFDSLDF  
DTLLKQAQKNLRR  
(1193-1542 aa encoded by BC063806)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.