

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SPDYA Fusion Protein



## Basic Information

**Catalog Number:**

Ag12041

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MRHNQMCETPPTVTVVVKSGSNRSHQPKKIPITLGRP  
ICKDNWQAF EKNTHTNNKSKRPKGPCLVIQRQDMTA  
FFKLFDDDLIQDFLWMDCCCKIADKYLAMTFVYFKR  
AKFTISEHTRINFFIALYLANTVEEDEEETKYEIFPWAL  
GKNWRKLFPNFLKLRDQLWDRIDYRAIVSRRCCCEVM  
AIAPTHYIWQRERSVHHSGAVRNYNRDEVQLPRGPS  
ATPVDCSLCGKKRRYVRLGLSSSSLSHTAGVTEKHS  
QDSYNSLSMDIIGDPSQAYTGSEVVNDHQSNKGGKIT  
NFLKKDKSMEWFTGSEE

(1-313 aa encoded by BC093005)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.