

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human MED27 protein



## Basic Information

**Catalog Number:**

Ag12046

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MADVINVSVNLEAFSQAISAIQALRSSVSRVFDCLKD  
GMRNKETLEGREKAFIAHFQDNLHSVNRDLNELERLS  
NLVGKPSENHPLHNSGLSLDPVQDKTPLYSQLQAYK  
WSNKLQYHAGLASGLLNQQLKRSANQMGVSAKRRP  
KAQPTTLVLPQQYVDDVISRIDRMFPEMSIHLSRPNGT  
SAMPLLVTLGKVLKIVVMRSLFIDRTIVKGYNENVYTE  
DGKLDIWSKSNYQVFQKVTDHATTALLHYQLPQMPD  
VVVRSFMTWLRSYIKLFQAPCQRCGKFLQDGLPPTW  
RDFRTLEAFHDTCRQ

(1-311 aa encoded by BC064608)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.