

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ADAL protein



## Basic Information

Catalog Number:

Ag12114

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MI E A E E Q Q P C K T D F Y S E L P K V E L H A H L N G S I S S H T M K  
K L J A Q K P D L K I H D Q M T V I D K G K R T L E E C F Q M F Q T I H  
Q L T S S P E D I L M V T K D V I K E F A D D G V K Y L E R S T P R R E N  
A T G M T K K T Y V E S I L E G I K Q S K Q E N L I D V R Y L I A V D R R  
G G P L V A K E T V K L A E E F F L S T E G T V L G L D L S G D P T V G Q  
A K D F L E P L L E A K I A G L K L A L H L S E I P N Q K K E T Q I L L D L L P  
D R I G H G T F L N S G E G G S L D L V D F V R Q H R I P L G K A W S F R  
S S R

(1-267 aa encoded by BC075857)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.