

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KCNAB2 protein



## Basic Information

Catalog Number:

Ag12120

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MYPESTTGSPARLSLRQTGSPGMIYSTRYGSPKRQLQF  
YRNLGKSGLRVSLGLGTWVTFGGQITDEMAEQMLTL  
AYDNGINLFDTAEVYAAAGKAEVVLGNIIKKKGWRRSS  
LVITTKIFWGGKAETERGLSRKHIIIEGLKASLERLQLEYV  
DVVFANRPDPNTPMEETVRAMTHVINQGMAMYWG  
TSRWSSMEIMEAYSVARQFNLTPPICEQAEYHMFQR  
EKVEVQLPELFHKIGVGAMTWSPLACGIVSGKYDSGI  
PPYSRASLKGQWLKDKILSEEGRRQQAKLKELQAIIE  
RLGCTLPQLAIAWCLRNEGVSLLGASNADQLMENE  
GAIQVRVRGPAGQRAHPSPSPVQCILPGSSCVPGSV  
LGTQDAPVNHQSCAPGELAFQEQEQT  
(1-395 aa encoded by BC110351)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.