

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# KCNAB2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag12120

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MYPESTTGSPARLSLRQTGSPGMIYSTRYGSPKRQLQF  
YRNLGKSGLRVSLGLGTWVTFGGQITDEMAEQMLTL  
AYDNGINLFDTAEVYAAAGKAEVVLGNIIKKKGWRRSS  
LVITTKIFWGGKAETERGLSRKHII EGLKASLERLQLEYV  
DVVFANRPDPNTPMEETVRAMTHVINQGMAMYWG  
TSRWSSMEIMEAYSVARQFNLTPPICEQAEYHMFQR  
EKVEVQLPELFHKIGVGAMTWSPLACGIVSGKYDSGI  
PPYSRASLKGQWLKDKILSEEGRRQQA KLKELQAI AE  
RLGCTLPQLAIAWCLRNEGVS SVLLGASNADQLM ENI  
GAIQVRVRGPAGQRAHPSPSPVQCILPGSSCVPGSV  
LGTQDAPVNHQSCAPGELAFQEQE T  
(1-395 aa encoded by BC110351)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.