

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human CD68 protein



Basic Information

Catalog Number:
Ag1215

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MRLAVLFSGALLGLLAAQGTGNDCPHKISATLLPSFT
VTPTVTESTGTTSHRTTKSHKTTTTHRTTTTGTTSHGPT
TATHNPTTSHGNVTVHPTSNSSTATSQGPSTATHSP
ATTSHGNATVHPTSNSSTATSPGFTSSAHPEPPPPSPS
PSPTSKETIGDYTWNGSQPCVHLQAQIQIRVMYTT
QGGGEAWGISVLNPNKTKVQGSCEGAHPHLLSFY
GHLSFGFMQDLQKVVYLSYMAVEYNVSPHAAQW
TFSAQNASLRDLQAPLQQSFSCSNSSII LSPAVHLDLL
SLRLQAAQLPHTGVFGQSFSCPSDRSILLPLIIGLILLGL
LALVLI AF CI IRRRPSAYQAL
(1-354 aa encoded by BC015557)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.