

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TIGD3 Fusion Protein



## Basic Information

**Catalog Number:**

Ag12253

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MELSSKKKLHALSLAEKIQVLELLDESKMSQSEVARRF  
QVSQPQISRI CKNKEKLLADWCSTANRERKRKRESK  
YSGIDEALLCWYHIARAKAWDVTGPMLLHKAKELADI  
MGQDFVPSIGWLVRWKRNNVGF GARHVLAPSFPE  
PPPPGLTSAQLPLSLKDFSPEDVFGCAELPLLYRAVP  
GSFGACDQVQVLLCANSRGTEKRRVLLGGLQAAPRC  
FFGIRSEALPASYHPDLGIPWLEWLAQFDRDMGQQG  
RQVALLLAARVVEELAGLPLGYHVKLLPLAASSTTPPLP  
SSVVRAFKAHYRHRLLGKLAAIQSERDGTSLAEAGAGI  
TVLDALHVASAAWAKVPPQLIFSSFIQEG LAPGKTPP  
SSHKTSEMPVPGGLSLEEF5RFVDLEGEPRSGVCKE  
EIGTEDEKGDREGAFELPTKADALRALGTLRRWFECN  
STSPELFEKFDCEEEVERLCCL  
(1-471aa encoded by BC074862)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.