

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human IGSF1 protein



## Basic Information

**Catalog Number:**

Ag12286

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

VTLWCRGPVHGVGYLHKEGEATSMQLWGSTSNDG  
AFPITNISGTSMGRYS CCYHPDWTSSIKIQPSNTLELL  
VTGLLPKPSLLAQPGPMVAPGENMTLQCQGELPDSTF  
VLLKEGAQEPLEQQRPSGYRADFWMPAVRGEDSGIY  
SCVYYLDSTPFAASNHSLSLEIWWTDKPPKPSLSAWP  
STMFKLGKDITLQCRGPLPGVEFVLEHDGEEAPQQFS  
EDGDFVINNVEGKIGNYSCSYRLQAYPDIWSEPSDP  
LELVGAAGPVAQECTVGNIVRSSLIVVVVALGVVLA  
IEWKKWPRLRTRGSETDGRDQTIALEECNQEGEPGTP  
ANSPSSTSQRISVELVPI

(987-1336aa encoded by BC063884)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.