

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human LYN protein



## Basic Information

**Catalog Number:**

Ag12310

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

KDITRKDAERQLLAPGNSAGAFLIRESLTKGSFSLSVR  
DFDPVHGDVIKHYKIRSLDNGGYISPRITFPICISDMIK  
HYQKQADGLCRRLEKACISPKPQKPWDKDAWEIPRE  
SIKLVKRLGAGQFGEVWMGYNNSTKVAVKTLKPGT  
MSVQAFLEENLMKTLQHDKLVRLYAVVTREPIYIIT  
EYMAKGSLLDFLKSDEGGKVLPLKLI D FSAQIAEGMAY  
IERKNIYHRDLRAANVLVSELMCKIADFGLARVIEDN  
EYTAREGAKFPIKWTAPEAIFGCFTIKSDVWSFGILL  
YEIVTYGKIPYPGRTNADVMTALSQGYRMPRVENC PD  
ELYDIMKMCWKEAEERPTFDYLSVLD D F Y T A T E G Q  
YQQQP

(132-512aa encoded by BC075002)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.