

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human INPP5F protein



Basic Information

Catalog Number:

Ag12323

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

VKQTKSNVNI~~GNLRKLG~~NFTKPEMKV~~NFLKPNL~~KVNL
WKS~~DSSLET~~MENTGVMDKVQAESD~~GMSSD~~NDSYH
S~~DEFLT~~NSK~~SDEDR~~QLANSLESVGPIDYVLPSCGIIAS
APRLGSR~~SQSL~~SST~~DSSV~~HAPSEITVAHGSLGK~~GQE~~
SPLK~~KSPS~~AGDVHILTGF~~AKPMD~~IYCHRFVQDAQNKV
THLSE~~TRSV~~SQQASQERNQMTNQV~~DETQ~~SESTEQ
TPSR~~PSQL~~DVSLSATG~~PQFL~~SVEPAHSVA~~SQK~~TP~~SA~~
SSM~~LELET~~GLHVTPSP~~ESSSS~~RAVSPFAKIRSSMVQV
ASITQAGLTHGINFAVSKVQKSPPEPEIINQVQ~~QNEL~~
KKMFIQCQTRIIQI

(C-term-345aa encoded by BC111493)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.