

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

INPP5F Fusion Protein



Basic Information

Catalog Number:

Ag12323

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

VKQTKSNVNI~~GNLRKLG~~NFTKPEMKVNF~~LKPNL~~KVNL
WKS~~DSSLET~~MENTGVMDKVAESD~~GDMS~~SDNDSYH
SDEF~~L~~TNSK~~S~~DE~~DRQLANS~~LESVGPIDYVLPSCGIIAS
APRLGSR~~SQ~~SL~~S~~ST~~DSSV~~HAPSEITVAHG~~S~~GLGK~~G~~QE
SPLK~~K~~SPSAGDVHILTGF~~AK~~PMDIYCH~~R~~FVQDAQ~~N~~KV
THLSE~~T~~RSV~~S~~Q~~Q~~AS~~Q~~ERN~~Q~~MTN~~Q~~V~~S~~DET~~Q~~SE~~T~~EQ
TPSR~~P~~SQ~~L~~DV~~S~~LSATG~~P~~Q~~F~~LS~~V~~EP~~A~~HS~~V~~AS~~Q~~K~~T~~PT~~S~~A
SS~~M~~LE~~L~~TGL~~H~~VT~~P~~SP~~S~~ESS~~S~~RA~~V~~SP~~F~~AK~~I~~R~~S~~SM~~V~~Q~~V~~
ASITQAGL~~T~~H~~G~~IN~~F~~AV~~S~~K~~V~~Q~~K~~SP~~P~~E~~P~~E~~I~~I~~N~~Q~~V~~Q~~N~~EL
KKMFI~~Q~~C~~Q~~TRII~~Q~~I

(C-term-345aa encoded by BC111493)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.