

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human DKC1 protein



## Basic Information

**Catalog Number:**  
Ag1239

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
VRLHNAIEGGTQLSRALETLTGALFQRPLIAAVKRQL  
RVRTIYESKMI EYDPERRLGIFWVSC EAGTYIRTL CVHL  
GLLGVGGQMQLRRVRS GVMSEKDHMVTMHDVLD  
AQWLYDNHKDES YLRRFVYPLEKLLTSHKRLVMKDSA  
VNAICYGAKIMLPGLRYEDGIEVNQEIVVITTKGEAI  
CMAIALMTTAVISTCDHGIVAKIKRVI MERDTPRKY  
GLGPKASQKMLMIKQGLLDKHKPTDSTPATWKQEVV  
DYSES AKKEVVAEVVKA PQVVAEAAKTAKRKRSESE  
ESDETPPAAPQLIKKEKKKSKDKKAKAGLES GAEPGD  
GDSDTTKKKKKKKKAKEVELVSE  
(C-term-350aa encoded by BC009928)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.