

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human RFTN2 protein



## Basic Information

**Catalog Number:**

Ag12397

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

CGLRKLEDDSSPGKIFSTLKRQVETKTEFAIEYVLL  
DFTLQASSNPEVIKINSILDIVTKVENYYLKGIVGAIH  
PVIQPVGQRKHLPAASYLYRVLLRLKLSPKNSAAPSGQ  
RRPRLVIEECPLTSEAQTNDAAKELIEKINVAARKGMK  
FVGFISQHYSPSKFCNGTNHDGDIESMLHVRHGSDE  
NCRSWNEGLSGQSSGESGIEELHHESGQYQMEQNG  
SPTSSKSRKGEASDNKLYTFNAFDDDDSTSWAYQEGI  
LSMKVTRKGSVISTLDADWLELTFYKQGLSLIDSFV  
FWETSKGEHLPKSLEGFFIYEEEGSGVPGSSRKGNDAI  
VVEQWTVIEGCEIKTDYGPLLHTLAIEFGWLLTSVLP  
VLRHDSGNLTKQIVFLQRPVMMNSAAQTPDKKAS  
RHIKGEDKNKATRSRIGLDTTSSQPAESRHLPEECRLSP  
SRECVTKGRLAQHNSFSGFSSSDNVLRELDGQFD  
QEDGVTQVTCM

(3-501 aa encoded by BC111069)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.