

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SFRS11 protein



## Basic Information

**Catalog Number:**  
Ag12457

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
MSNTTVVPSTAGPSPGGPGGGGGGGGGGGTEVI  
QVTNVSPSASSEQMRTLFGFLGKIDELRLFPPDDSP  
VSSRVCFVKFHDPSAVVAQHLTNTVFVDRALIVVP  
YAEGVIPDEAKALLAPANAVAGLLPGGGLLPTNPL  
TQIGAVPLAALGAPTLDPALALGLPGANLNSQSLAA  
DQLLKLMSVDPKLNHVAAGLVSPSLKSDTSSKEIEEA  
MKRVREAQSLISAAIEPDKKEKRRHSRSRSRRRTP  
SSSRHRRSRSRRRSHSKRSRRRSKSPRRRSHSRER  
GRRSRSTKTRDKKEDKPKRSTPPKSYSTARRSRA  
SRERRRRSRSGTRSPKPRSPKRKLSRSPRRHKKEK  
KKDKDKERSRDERERSTSKKKKSKDKDRERKSESDK  
DVKVTRDYDEEEQGYDSEKKEKKEKPIETGSPKTKEC  
SVEKGTGDSLRESKVNDDHHEEDMDMSD  
(1-483aa encoded by BC040436)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.