

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SFRS11 Fusion Protein



Basic Information

Catalog Number:
Ag12457

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MSNTTVVPSTAGPGSPGGPGGGGGGGGGGGTEVI
QVTNVSPASSEQMRTLFGFLGKIDELRFPDDSP
VSSRVCFVKFHPDSAVVAQHLTNTVFVDRALIVVP
YAEGVIPDEAKALSLAPANAVAGLLPGGGLLPTNP
TQIGAVPLAALGAPLDPALALGLPGANLNSQSLAA
DQLLKLMSVDPKLNHVAAGLVSPSLKSDTSSKEIEEA
MKRVREAQSLISAAIEPDKKEKRRHSRSRSRRRTP
SSSRHRRSRSRRRSHSKRSRRRSKSPRRRSHSRER
GRRSRSTSKTRDKKEDKSKRSTPPKSYSTARRS
SRERRRRSRSGTRSPKPRSPKRKLSRSPRRHKKEK
KKDKKERSRDERERSTSKKKKSKDKEKDRERKES
DVKVTDRDYDEEEQGYDSEKKEKKEKPIETGSPK
SVEKGTGDSLRESKVNDDHHEEDMDMSD
(1-483aa encoded by BC040436)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.