

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HS3ST5 Fusion Protein



## Basic Information

**Catalog Number:**

Ag12532

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MLFKQQA~~WLRQ~~KLLVLGSLAVGSLLYLVARVGS~~LDRL~~  
QPICPIEGR~~LGGARTQAEF~~PLRALQFKRGLLHEFRKGN  
ASKEQVRLHDLVQQLPKAIIIGVRKGGTRALLEMLNLH  
PAVVKASQEI~~HFFD~~NDENYKGI~~EWYR~~KMPFSYPO  
QITIEKSPAYFITEEVPERIYKMNSSIKLLIIVREPTTRAI  
SDYTQVLEGKERKNKTYKFEKLAIDPNTCEVNTKYKA  
VRTSIYTKHLERWLKYFPIEQFHVVVDGDRLITEPLPELQ  
LVEKFLNLPPRISQYNLYFNATRGFYCLRFNIIFNKCLA  
GSKGRIHPEVDPSVITKLRKFFHPFNQKFYQITGRTLN  
WP

(1-346aa encoded by BC093911)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.