

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ICA1L protein



Basic Information

Catalog Number:

Ag12570

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

EVATFSQRAVSDTLMTINRMEQARTEYRGALLWMKD
VSQELDPDTLKQMEKFRKVMQVRNSKASFDKLMKD
VCQKVDLLGASRCNMLSHSLTTYQRTLLGFWKKTAR
MMSQIHEACIGFHPYDFVALKQLQDTPSKI SEDNKDE
QIGGFLTEQLNKLVLSD E E A S F E S E Q A N K D H N E K H S Q
MREFGAPQFSNSENVAKDLPVDSLEGEDFEKEFSFLN
NLLSSGSSSTSEFTQECQTAFGSPSASLTSQEPSMGS
EPLAHSSRF LPSQLFDLGFHVAGAFNNWVSQEESELC
LSHTDNQPVPSQSPKKLTRSPNNGNQDMSAWFNLF A
DL D P L S N P D A I G H S D D E L L N A

(133-482 aa encoded by BC063876)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.