

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human HS3ST5 protein



## Basic Information

Catalog Number:

Ag12621

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MLFKQQA~~WLR~~QKLLVLGSLAVGSLLYLVARVGS~~LDRL~~  
QPICPIEGR~~LG~~GARTQAEFPLRALQFKRGLLHEFRKGN  
ASKEQVRLHDLVQQLPKAIIIGVRKGGTRALLEMLNLH  
PAVVKASQEI~~HF~~FDNDENYKGI~~EWYR~~KMPFSYPO  
QITIEKSPAYFITEEVPERIYKMNSSIKLLIIVREPTTRAI  
SDYTQVLEGKERKNKTYKFEKLAIDPNTCEVNTKYKA  
VRTSIYTKHLERWLKYFPIEQFHVVVDGDRLITEPLPELQ  
LVEKFLNLPPRISQYNLYFNATRGFYCLRFNIIIFNKCLA  
GSKGRIHPEVDPSVITKLRKFFHPFNQKFYQITGRTLN  
WP

(1-346aa encoded by BC093911)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.