

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human tubulin- Alpha protein



## Basic Information

Catalog Number:

Ag12654

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MRECISVHVGGQAGVQIGNACWELFCLEHGIQADGT  
FDAQASKINDDDSFSTFFSETGNGKHVPRAVMIDLEP  
TVVDEVVRAGTYRQLFHPEQLITGKEDAANNYARGHY  
TVGKESIDLVLDRIRKLTACSGLQGFLIFHSFGGGTG  
SGFTSLLMERLSLDYGKSKLEFAIYPAPQVSTAVVEP  
YNSILTTHTTLEHSDCAFMDNEAIYDICRRNLDIERPT  
YTNLNLISQIVSSITASLRFDGALNVDLTFQTNLVPY  
PRIHFPLVTYAPIISAEKAYHEQLSVAEITSSCFEPNSQ  
MVKCDPRHGKYMACCMLYRGDVVPKDVNVAIAAIKT  
KRTIQFVDWCPTGFKVGINYQPPTVVPGGDLAKVQR  
AVCMLSNTTAIAEAWARLDHKFDLMYAKRAVHWY  
VGEGMEEGEFSEAREDLAALEKDYEEVGTDSFEEENE  
GEEF

(1-449aa encoded by BC074827)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.